Dentsply Sirona Universal Stain and Glaze are indicated for staining and glazing the following restorations:

- > Celtra Press zirconia reinforced lithium silicate (ZLS)
- > Celtra DUO zirconia reinforced lithium silicate (ZLS)
- > Cercon ht zirconia
- > Ceramco PFZ porcelain for zirconia
- > Cercon Ceram Kiss porcelain
- > Ceramco3 porcelain
- > Ceramco iC porcelain
- > Duceram Kiss porcelain

Contraindications

Only the indications listed above are suitable.

Warnings/Adverse Reactions

Consult the Dentsply Sirona Universal Stains and Glaze Safety Data Sheets (SDS) for safe handling and usage. Use only according to the indications listed above.

Undesired side effects of these medical products are extremely rare with proper processing and use. In case of skin sensitization or rash, discontinue use and seek medical attention. Immunoreactions, such as allergies and/or local irritations (taste or oral mucosa), cannot completely be excluded on principle.

For patients with hypersensitivity to any of the ingredients, this medical device must not be used or only used under strict supervision of a treating physician/dentist. Similarly, known cross-reactions or interactions of this medical product with other materials existing in the mouth must be considered by the physician/dentist.

Precautions

- > Wear suitable protective eyewear, clothing and gloves
- > Do not inhale abrasive dusts
- > Do not ingest any of the materials
- > Frameworks that show evidence of cracking or surface voids must be discarded and not processed further

Technical data

- > Coefficient of Thermal Expansion (CTE): $8.0 \pm 0.5 \times 10^{-6} \text{K}^{-1}$ (25-400°C)
- > Device classification: Dental ceramics, type I, class 1b per ISO 6872: 2015
- > Chemical solubility: $10 \mu\text{g}/\text{cm}^2$

Storage Conditions

Avoid direct exposure to sunlight and high heat, protect from moisture, when not in use keep lid tightly closed.

Symbols on product labels

	Product code or reorder number
	Batch or Lot number
	Expiration date
	Keep dry
	Date of manufacture
	Consult instructions for use

Use Instructions:

Framework Preparation:

Follow the current corresponding Directions for Use for preparing the substructure material.

Directions for Use:

- > Celtra Press zirconia reinforced lithium silicate - 50539994
- > Celtra DUO zirconia reinforced lithium silicate - 50539985
- > Cercon ht zirconia - 090961
- > Ceramco PFZ porcelain for zirconia - 090862
- > Cercon Ceram Kiss porcelain - 28078
- > Ceramco3 porcelain - 301700
- > Ceramco iC porcelain - 745-215
- > Duceram Kiss porcelain - 22147

Stain & glaze application:

Note: If the stain or glaze inside the jar has separated, mix thoroughly with a glass, plastic or zirconia spatula. Withdraw desired amount of Universal Stain or Glaze from the jar and place it on the mixing palette. If a thinner consistency is desired, dilute the material with the Dentsply Sirona Stain and Glaze Liquid (REF 601315/15mL, 601350/50mL). (*See Figure 1*)

Apply a thin layer of the stain or glaze material all over the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin. Too thin of a layer will result in a more matte finish. Too thick of a layer can result in puddling and/or pitting of the material. (*See Figure 2*)

If a more intensive shade effect is desired, it is achieved by several staining applications and firings. However, the application of too many stain layers can result in an unnatural shade effect.


Figure 1

Figure 2

Stain & glaze firing:

Please use the suggested furnace cycle/program shown below:

	Drying (min)	Pre-heat (min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air (min)	Cool Time (min)
Celtra Press (full contour)	2	2	400°C	55°C/min	760°C	0	0	0	2:00	5
Celtra Press (cutback or 2nd firing of full contour)	2	2	400°C	55°C/min	750°C	0	0	0	2:00	5
Celtra DUO	2	2	500°C	55°C/min	820°C	0	0	0	1:30	3
Cercon ht	3	3	450°C	55°C/min	820°C	0	0	0	1:30	6
Ceramco PFZ	5	5	450°C	60°C/min	850°C	0	0	0	0:30	6
Cercon Ceram Kiss	2	3	450°C	55°C/min	800°C	0	0	0	1:00	6
Ceramco3	3	3	650°C	70°C/min	935°C	0	0	0	0:30	0
Ceramco iC	3	3	500°C	35°C/min	750°C	500°C	750°C	1	1:30	0
Duceram Kiss	3	1	575°C	55°C/min	890°C	0	0	0	1:00	0

* Drying time depends on the type of furnace

* Additional Stain/glaze and characterization firing cycles can be conducted with the same firing parameters. Do not go below or above the temperature range 740°C to 950°C.

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DENTSPLY International Inc.
York, PA 17401 U.S.A.