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Ducera® Lay SUPERFIT

Ducera® Lay SUPERFIT

1. Prepare the master or the working models.
2. Block out undercuts.
3. Apply an even layer of Ducera® Lay distance laquer (spacer) to the die as space for bonder or cement, but leave out an area of approximately 1 mm around the margin.

As soon as the surface is dry you can start the duplication.

4. Place the master model in the duplication flask.
5. We recommend the use of addition cross linked silicon as duplicating material.
6. The correct mixing ratio is essential for a good result:

Mixing Ratio: 10 parts of refractory powder : 3.5 parts of pure mixing liquid (Important: SUPERFIT mixing liquid is a ready-to-use liquid!

Do not add any water or other liquids!). In order to achieve steady results, the exact mixing ratio of SUPERFIT refractory powder and SUPERFIT mixing liquid should always be determined with a scale (gold scale with 0.01 g accuracy).

- To achieve a homogenous mixture it is necessary to mix a minimum of 30 g refractory powder and 10.5 g mixing liquid.
- Adding a larger amount of mixing liquid improves the flow of the mixture, however, at the same time decreases the final strength and alters the expansion characteristics. Please note that the use of a wet mixing bowl can also affect the expansion.

7. Premix the Ducera® Lay SUPERFIT by hand for a short time. Be careful not to create any bubbles.
- 8 a. Then mix it in a vacuum mixer for about 45 sec. (Important: Using less than 30 g SUPERFIT powder results in an inhomogenous mixture).
- 8 b. If a small vacuum mixer is not available, mix intensively by hand and then place the mixture into a vacuum container for 30 sec. to degas it.
9. Place the impression form on a vibrator to pore the refractory die into the form. Avoid bubbles!
10. Remove the die after 45 mins (refer to the instructions of the duplicating material manufacturer).

SINTER FIRING: Ducera® Lay SUPERFIT dies are sintered as follows:

1. Pre-heat at 575° C for about 5 mins.
2. Raise the temperature by 55°C/min. up to 1100°C without vacuum.
3. Hold the final temperature for about 5 mins.

When sintering more than 5 dies at the same time, increase the holding time to 7 mins.

The same applies for sintering in the Systemat M/D.

After sinter firing, the refractory dies should be white. A discoloration into grey is a sign that:

- the sinter temperature is too low
- the holding time is too short

Sinter firing compensates for the setting expansion.

If the above mentioned instructions are followed correctly, the Ducera® Lay SUPERFIT die and the master model will be of exactly the same size.

Firing Chart for Ducera® Lay Connector Paste:

Predrying:	1 min.	Heat rate:	55°C with vacuum
Preheating:	3 mins.	Final temperature:	970°C
Preheating temperature:	620°C	Holding time:	1 min.

Firing Chart for Duceram® metal ceramic on the refractory die:

Dentine firing: Use a final temperature of about 30°C above the regular metal ceramic's firing temperature.

Important: We recommend max. heat rate (approximately 130°C/min.) and fast cooling for all firings of Duceram® porcelain to refractory material.

Please note: In case refractory dies are not sintered immediately after removal from the duplication form, a velvet-like film can develop on their surface. This process is favored by high humidity. The film does not influence the surface characteristics and the fitting qualities of Ducera® Lay SUPERFIT. The velvet like film does not develop when the die is sintered immediately, respectively disappears after sintering.