



Transfer table

Interaction → Kiss



REV 2014-10

GB

Interaction	Kiss
 interaction®	
Incisals	
incisal 57	S 1
incisal 58	S 1 + S 2 (mix 1:1)
incisal 59	S 2
incisal 60	S 3
x-tra incisal clear	TC
x-tra incisal blue	Ocean (<i>use in diluted/attenuated form</i>)
x-tra incisal grey	Fog (<i>use in diluted/attenuated form</i>)
x-tra incisal orange	Sunrise
x-tra incisal red	TR
x-tra incisal white	WS
x-tra incisal medium	Clear + Sky (mix 1:1)
x-tra incisal bright	Flu 1 + Clear (mix 1:1)
x-tra incisal dark	OS 2
Chroma-supporting additional masses	
action-i dentin 1A2	PC 2
action-i dentin 3A4	PC 4
action-i dentin 1B2	P 1 + P 3 (mix 1:1)
action-i dentin 3B4	P 1 + P 4 (mix 1:1)
action-i dentin 1C2	P 4 + P 6 (mix 1:1)
action-i dentin 3C4	PC 6
Gingival masses (gums)	
gum 213 / light	Gum 1
gum 212 / dark	Gum 2
Paste opaque Gum	Paste opaque Gum